



# Contactless infrared forehead thermometer

Medical grade infrared frontal body measuring instrument

SUPER INSTRUMENT LTD

Accuracy. Reliability. Speciality





ACCURATE TEM  
MEASUREMENT

精准测温

进口优  
精准



## OUTLINE

01

Product Background

02

Product Features

03

Product Technical  
Parameters

04

Product Instructions



# Product Background

## Product Background



Infrared thermometer is a thermometer based on the principle of radiation to measure human body temperature. The infrared sensor used only absorbs the infrared radiation radiated from the human body without emitting any rays to the outside world and senses the human body temperature via the contactless method.



### Working Principle

The human body's infrared thermal radiation is focused on the detector, which converts the radiated power into an electrical signal. This signal can be displayed in units of temperature after being compensated by the ambient temperature.



### Advantages

High accuracy Real-time Fast





## Product Background



### Infrared thermometer

The infrared thermometer inputs the error calculation formula of the temperature difference into the central micro-processing system, detects the change of the ambient temperature, and corrects its reading value at any time, to ensure that the infrared thermometer measurement is more accurate. Infrared thermometers are mainly divided into contact and contactless types.



### Categories of infrared thermometer

**Contact infrared thermometer:** Taking the ear thermometer as an example, it is mainly composed of shell, temperature probe, temperature sensor, PCB, LCD, buzzer and battery

**Non-contact infrared thermometer:** Taking the frontal temperature gun as an example, it is mainly composed of infrared probe, main circuit board, LCD display and shell.





## Product Features



Accuracy. Reliability. Speciality

- High accuracy
- Fast
- Safe
- Non-contact

## Product Features



### Forehead thermometer

Based on the principle of infrared radiation, the non-contact design is relatively safe, accurate and thus suitable for the elderly, children, hospital and family.

## Forehead thermometer Introduction



Temperature  
32°C ~ 42.5°C

Net Weight  
136g (Without battery)

Size  
L100xW46xH160m  
m

Accuracy. Reliability. Speciality



## Forehead thermometer Introduction



"+" "-" switch

Page-turning Button

Review Button



**Remind function**

**32 sets of memory data**

Conveniently view  
comparative data and  
keep an eye on body  
temperature changes

## Forehead thermometer Features



### High accuracy

High quality probe with high-precision infrared sensor, Automatic measurement, accurate calculation, and real body temperature reflection



### Fast

The measurement can be easily completed in 0.5-1 second



### Safe

The instrument only absorbs the weak infrared measurements emitted by the human body, and the instrument itself does not emit any rays, which is safe and healthy.



### Non-contact

The non-contact design to eliminate cross-infection.

High accuracy



# ACCURATE TEMPERATURE MEASUREMENT

High-precision infrared sensing, automatic measurement, accurate calculation, and real body temperature reflection



红外额式体温计

型号: MD-33520 设备编号:  
测温范围: 32°C~42.5°C 测量误差:  $\pm 0.3^{\circ}\text{C}$   
工作电压: d.c. 9V(6F22)  
生产日期: 2020年3月3日  
使用期限: 3年



其他内容详见说明书



SECOND-SPEED TEMPERATURE  
MEASUREMENT REQUIRES NO EXTRA TIME.

**The measurement can be easily completed  
in 0.5-1 second**



Accuracy. Reliability. Speciality





## NO RADIATION SAFETY PROTECTION

**The instrument only absorbs the weak infrared measurements emitted by the human body, and the instrument itself does not emit any rays, which is safe and healthy.**





## »»» Contactless design



NON-  
CONTACT


Non-contact design to  
eliminate cross-infection,  
which is safer.

Accuracy. Reliability. Speciality



## Product Technical Parameters

## Product Technical Parameters



### 产品参数 Technical Parameters

Display resolution 0.1°C

Test error  $\pm 0.3^{\circ}\text{C}$

Measurement distance: 5-15cm

Operating environmental temperature: 10-40°C

Operating environment humidity:  $\leq 80\%$

Power: DC9V battery

Auto shut off: 10 seconds

Size 100x46x160mm(Lx WxH)

Weight 125g(Without battery)

Accuracy. Reliability. Speciality



### Comparison of different temperature measurement methods

The body temperature measured by different temperature measurement methods will be different from the specific temperature. The specific temperature differences are shown in the table below.

Forehead Temperature 35.8°C ~ 37.8°C

Cochlea Temperature 35.8°C ~ 38°C

Oral Temperature 35.5°C ~ 37.5°C

Oxter Temperature 34.7°C ~ 37.3°C

Anus Temperature 36.6°C ~ 38°C

## Product Technical Parameters



### A professional temperature measurement instrument

01

#### Temperature within 32°C ~ 42.5°C

The product displays temperature when the temperature is measured within 32°C ~ 42.5°C

02

#### Temperature > 42.5°C

The product will display HI when the temperature excess 42.5°C

03

#### Temperature < 32°C

The product will display LO when the temperature is lower than 32°C







## Product Instructions





## FRONTAL TEMPERATURE GUN

24 小时给宝宝测温  
给他最好的呵护



### 医用级 专业红外额温枪

A professional body temperature measuring instrument,  
always watching body temperature changes

## Product quick use guide

- ◆ The thermometer is in temperature measurement mode. Press and hold the power button to measure continuously. If you do not need a long measurement, please release it after pressing the power-on button.
- ◆ The counter "LOG" is behind the lower left corner of the screen. The count value is cycled from 1 to 32. The counter will increase by 1 each time a valid measurement is performed, and then return to 1 after 32. If the measurement result is "LO" or "HI", it is called invalid measurement, the measurement result is not entered into the memory, and the counter is not incremented.
- ◆ Historical data view. After each test, the instrument will automatically record the test data, up to 32 historical measurement data. To view the data, press "M" key and then press "+" or "-". Press "-" to display the last temperature measured. Press "+" is the farthest record. Press "+" or "-" again to loop.

## Button Instructions



01

Single-press "M" key for historical data query, press "+" to view the next group of temperature readings, and press "-" to view the previous group of temperature readings. A total of 32 sets of historical data can be viewed.

02

If the measurement result is "LO" or "HI", it is called invalid measurement, the measurement result is not entered into the memory, and the counter is not incremented.

03

Press the backlight button to turn the backlight on and off

04

Press the horn button to turn sound on and off



## Qualification



## Authorized certification Quality assurance

The product has national medical registration  
license and passed the CE/RoHS certification

Accuracy. Reliability. Speciality

# THANKS

For Your Attention



CORETERM

Camí Vell de Russafa, 713 Pol. Ind. Catarroja 46470 CATARROJA ( VALENCIA ) [coreterm@coreterm.es](mailto:coreterm@coreterm.es)